



Bariatric Products



Bariatric Products

To cater adequately for large and heavy patients, our design team has focused on larger, reinforced frameworks, electronics with increased weight capacity and specialised systems such as patient weigh scale and powered-drive options.

Specialised Bariatric Chairs and Products in Australia

Patients who are heavier than average can feel comfortable and stable with our selection of bariatric chairs in Australia. It's crucial for individuals in the medical industry, in particular, to attend to the individual needs of each potential patient. Pushing patients around in hospitals and other medical facilities no longer requires any human effort thanks to our bariatric chairs.



Contour Recline Barituff



Contour Energy



Contour Barituff



Contour Barituff-X



Orbit-Drive



Weigh Scales



Aquatuff Bari

Contour Recline Barituff

The Barituff is a 500kg SWL bariatric chair bed with an extra-wide top and deeper seat.



Key Benefits

Low & safe: The ultra-low height is combined with dropsides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.

Purpose built: Increased componentry strength in the core frame design operates with higher capacity columns to safely lift 500kg without overload.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Day Surgery, ED, Podiatry & General Ward

Why Buy?

Having equipment that is robust enough, that you can completely rely on, is the greatest challenge when caring for bariatric patients.

We have increased componentry strength in the core frame design which operates with higher capacity columns to safely lift 500kg without overload. A wider top with a deeper seat has been standardized to comfortably cater for large patients. Heavy-duty, dark coloured vinyl is used on the mattress cushions to combat wear and staining. This stretcher chair has been built for bariatric patients.

Contour Recline Barituff



Low and safe

Excellent 480mm low height with fold-under dropside.

The ultra-low height is combined with dropside that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



Easy egress

Forward tilt facilitates getting on and off.

For aged or disabled patients, a gentle incline and the long dropside armrests may assist in safer transfer.



Lies flat

The backrest and leg rest move independently or together from chair to bed.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors. The unique handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



Large top design

The Barituff has a 200mm wider x 140mm deeper seat than the standard Recline.

The seat size is increased to 850mm wide x 600mm deep, easily accommodating larger patients. Use our range of colour coding choices if required to assist in quick identification.



Clear handset

Large simple functions with colour assistance.

Position the patient quickly and easily with the touch of a button - either by independently moving sections or applying one touch contouring.

Specifications

Size: 980mm wide x 1270mm upright (1750mm lying flat) overall

Seat width: 850mm wide x 600mm deep seat

Mattress size: 110mm pressure care mattress cushion with heavy duty cover

Height Range: 480 - 880mm with lie-flat and forward/reverse tilt

SWL: 500kg

Castor size: Twin wheel 150mm castors

Package Includes: Standard IV, vertical oxy holder, drop down footrest & meal tray

Product Code

Contour Recline Barituff: MPC-C-REB

Contour Recline Barituff Package: MPC-C-REB-P

Contour Energy

A transport chair with ZERO effort, 4-way power drive reducing the risk of strain during patient transport.



Key Benefits

Ergonomics & intuitive: Drive effortlessly in any direction with only the slightest of touch, paired with perfect positioning of the drive controls for ease of use.

Easy access: The large footplates offer generous support to your patients and lift quickly for clear egress.

Powerful movement: A rotating electric drive wheel provides excellent traction, while center steering makes patient transport easier over longer distances.

Why Buy?

Take all the human effort out of pushing and maneuvering chairs with an intuitive and intelligent drive system designed for easy use by any operator. Natural movement controls the speed and direction of the chair with the ability to stop very quickly by pulling back. Help prevent the most common carer strain injuries by taking the back twist out of moving heavy weight forwards, backwards and sideways. Patients are positioned comfortably with a gentle recline and feet fully supported. We've taken care with every aspect of design to ensure comfort and efficiency for both the patient and carer. Customise for a wide range of healthcare areas with relevant accessories.

Contour Energy



Ergonomic

Drive effortlessly in any direction by only light touch.



Intuitive

Perfect positioning of the drive controls for the operator.



Support

The large foot plates offer generous support to the patient.



Easy access

The large foot plates lift quickly and easily for clear egress.



Full clearance

Complete frontal clearance by lifting and rotating foot plates outwards.



Removable armrests

The side armrests can be fully removed for multi direction egress.



Versatile

Use the storage shelf during transport or customise with accessories such as IV poles, oxy holders and seat belt.

Specifications

Size: 750mm W x 1300mm L x 1180mm H

Seat Size: 650mm W x 520mm D x 550mm H

SWL: 300kg

Castor Size: 150mm

Package Includes: standard IV & vertical oxy holder

Product Code

Contour Energy: MPC-C-ENT-EDU

Contour Energy Package: MPC-C-ENT-EDU-P

Contour Barituff

A 500kg SWL extra-wide electric stretcher designed for Bariatric patients.



Key Benefits

Purpose built: Increased componentry strength in the core frame design operates with higher capacity columns to safely lift 500kg without overload.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Day Surgery, ED, General Ward & Bariatric Care

Why Buy?

Having equipment that is robust enough, that you can completely rely on, is the greatest challenge when caring for Bariatric patients. The Contour Barituff has increased componentry strength in the core frame design and operates with higher capacity columns to safely lift 500kg without overload.

The stronger and wider top comfortably caters for large patients, with a choice of mattress styles available to suit the type of departmental use. Stretcher stability is assured through an increased base size. Powered drive and weigh scales genuinely reduce manual handling risks to the carer. This stretcher has been built for Bariatric patients.

Contour Barituff



Ultra-low height

Easier for patients to get on and off.
The 470mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.
Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 850mm (plus an extra 110mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

Specifications

Overall Size: 2135 x 1100mm

Patient Surface: 1900 x 980mm

Mattress: 110mm thick with 150mm thick option

Height Range: 470 - 850mm

SWL: 500kg

Castor Size: 150mm twin wheel castors

Package Includes: 110mm premia bari mattress, folding IV, vertical oxy, weigh scales, orbit drive, emergency CPR backrest release, 2nd handset & USB port

Product Code

Contour Barituff: MPT-C-BE

Contour Barituff Package: MPT-C-BE-P

Contour Barituff-X

A 500kg SWL extra-wide electric trauma stretcher designed for Bariatric patients.



Key Benefits

Purpose built: Increased componentry strength in the core frame design operates with higher capacity columns to safely lift 500kg without overload.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for ED / Imaging and Bariatric Care

Why Buy?

Having equipment that is robust enough, that you can completely rely on, is the greatest challenge when caring for Bariatric patients. The Contour Barituff-X has increased componentry strength in the core frame design and operates with higher capacity columns to safely lift 500kg without overload.

The stronger and wider top comfortably caters for large patients, with a choice of mattress styles available to suit the type of departmental use. Stretcher stability is assured through an increased base size. Powered drive and weigh scales genuinely reduce manual handling risks to the carer.

This stretcher has been built for Bariatric patients.

Contour Barituff-X



Ultra-low height

Easier for patients to get on and off.
The 600mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.
Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 980mm (plus an extra 110mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Cassette holder options

Holders can be produced with single or double sided handles and any depth or shaped rebates to suit all of the unique cassette sizes in use today

Specifications

Overall Size: 2135 x 1100mm

Patient Surface: 1900 x 980mm

Mattress: 110mm thick with 150mm thick option

Height Range: 600 - 980mm

SWL: 500kg

Castor Size: 150mm twin wheel castors

Package Includes: 110mm Premia Bariatric mattress, folding IV, weigh scales, vertical oxy, orbit drive, 2nd handset & USB port

Product Codes

Contour Barituff-X: MPT-C-BEX

Contour Barituff-X Package: MPT-C-BEX-P

X-Ray Cassette Holder: MPT-CA-XCH-D

Orbit-Drive

Taking patient transport to a new level in carer welfare with zero effort, 4-way power drive.



Key Benefits

- Forward, back and sideways powered drive
- Intuitive height adjustable controls with emergency stop
- Powerful electric drive centre steer wheel

Orbit-Drive



Adjustable controls

Lower the drive controls down below the mattress height for enhanced patient access and specialised accessory use.



Adaptable

Install the power drive system on many of the Contour stretcher models



Powerful drive

An electric driven wheel pressured onto the floor that provides excellent traction and center steering.



Sideways movement

The unique rotating design allows the stretcher to be driven in any direction - forwards, backwards and sideways!

Why Buy

Take all the human effort out of pushing and jostling stretchers with an intuitive and intelligent drive system designed for easy use by any operator. Natural movement controls the speed and direction of the stretcher with the ability to stop very quickly by pulling back.

Help prevent the most common carer strain injuries by taking the back twist out of moving a stretcher forwards, backwards and sideways. Importantly, we haven't overlooked safety features such as emergency stop and system lock out, which are standard. Add this system to our traditional vertical column stretcher designs such as Contour Barituff, Classic, Deluxe and Multi-X to improve both standard and Bariatric patient transport.

Product Code

- MPT-CA-EDU

Weigh Scales

Add a weigh scale system to Contour Stretchers for ease of patient measuring and tracking.

Key Benefits

200g accuracy Auto off with weight memory function Freeze weight to allow for additional accessories Excludes Contour Portare stretcher models

Product Code

- MPT-CA-WS



Easy to use

Aquatuff Bari

A robust, fully electric 300kg SWL shower trolley with extra wide base and top for Bariatric patients.



Key Benefits

Safe & secure: The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfers.

Hygiene conscious: Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.

Custom adjustments: Raise your patient to a comfort position with a 15 degree electric backrest option available.

Application

Designed for heavy use and demanding environments.

Why Buy?

Having equipment that is robust and completely reliable is of the utmost importance when caring for larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the lifetime of value that this bariatric shower trolley will give. The electric operation helps to safely lift 300kg without overload, whilst four twin wheels with brakes with wider base and top all increasing versatility in special needs.

Aquatuff Bari



Safe transfer

The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfer.



Easy drainage

Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.



Backrest option

Raise your patient to a comfort position with a 15 degree electric backrest.



Length adjustment

Lock the removable end rails (on each end of the top) into two positions, allowing 120mm extra space for longer patients.

Specifications

Size: 1195 W x 1950mm including bumpers

Liner Inner (min) Base Size: 950 W x 1700 L x 180mm Deep

Standard Hose Length: 1100mm, custom lengths available

SWL: 300kg

Height Range: 680 to 900mm

Castor Size: Quality 125mm twin wheel castors with total lock and steer

Electric Backrest Option: 15 degree

Construction: Powder coated stainless steel

Product Code

Aquatuff Bari; MST-BARI

Electric backrest option - MST-SA-EBB

Notes

Notes
