

# Contour Endosurge

An endoscopy imaging stretcher with a sliding top for quick and easy patient positioning.



## Key Benefits

**Smooth and steady:** An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor.

**Effortless control:** Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

**Supercharged:** Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

**Longest warranty:** Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

## Application

Ideal for Endoscopy / Colonoscopy & Pain Management

## Why Buy?

Take away the effort of re-positioning a heavy C-Arm X-Ray camera. For quick and easy movement of the patient, unlock the stretcher top and slide it a generous distance back or forwards.

Our aim is to provide more freedom and take strain off your backs. The cutaway base cover allows more clearance for the camera under the stretcher whilst bringing it closer to the floor for a lower working height.

v-01032023

# Contour Endosurge



## Slide to head

Unlock and move the top 140mm towards the head end from centre position.

Release the top lock and easily slide the stretcher top back or forwards to position your patient without having to significantly move the heavy C-arm imaging camera.



## Slide to foot

Unlock and move the top 260mm towards the foot end from the centre position.

Release the top lock and easily slide the stretcher top back or forwards to position your patient without having to significantly move the heavy C-arm imaging camera.



## Hydraulic option

Non-powered, manually operated stretcher option.

In some cases, such as small or remote facilities, or large high pressure departments, hydraulically operated stretchers may be more beneficial, with simple mechanical servicing and without the requirement of regular charging or power cords.

## Specifications

**Overall Size:** 2135 x 815mm (705mm narrow & 965mm wide options)

**Patient Surface:** 1900 x 680mm (570mm narrow & 830mm wide options)

**Electric Height Ranges:** 570 - 940mm

**Hydraulic Height Ranges:** 630 - 960mm

**SWL:** 300kg

**Castor Size:** 200mm castors and Ezi-Trak 5th wheel steering

**Sliding Top:** 400mm slide; 140mm to head and 260mm to foot from centre

**Electric Package Includes:** 100mm premia mattress, folding IV, vertical oxy, dropdown footboard, drop down push handles, emergency CPR backrest release, split dropsides, 2nd handset & USB port

**Hydraulic Package Includes:** 100mm premia mattress, folding IV, vertical oxy, dropdown footboard, drop down push handles & split dropsides

**Options:** Split dropside option (shown) and mid-position catch upgrades

## Product Code

Electric: MPT-C-EE

Electric Package: MPT-C-EE-P

Hydraulic: MPT-C-EH

Hydraulic Package: MPT-C-EH-P

v-01032023