

Contour Heli-Transfer

A stretcher designed to reduce helicopter patient transfer effort outdoors.



Key Benefits

Smooth and steady: An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for patient transfers

Why Buy?

Safely and easily move a helicopter stretcher across onto this transfer stretcher, a platform designed to minimise effort whilst transporting patients over rougher outdoor surfaces. After benefiting from the smooth, flat bridge that lines up with the helicopter door, rapidly fasten the helicopter stretcher securely down with the four quick pull straps.

Reduce the trauma to patients and babies in cribs and the workload for carers whilst transferring via the large diameter castors and very stable stretcher.

v-17012024

Contour Heli-Transfer



Hydraulic option

Non-powered, manually operated stretcher option.

In some cases, such as small or remote facilities, or large high pressure departments, hydraulically operated stretchers may be more beneficial, with simple mechanical servicing and without the requirement of regular charging or power cords.



Orbit Drive Upgrade

Help prevent strain injury with this unique forward, back and sideways powered drive system.



RACQ Life Flight

Built in backrest with a special thin high density mattress and split dropsides.



Westpac

6 riser dopside and rubberised outdoor castor wheels.

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 700mm (590mm narrow & 850mm wide options)

Electric Height Ranges: 550 - 920mm

Hydraulic Height Ranges: 530 - 860mm

SWL: 300kg

Castor Size: 200mm rubberised castors for outdoor shock absorption

Package Includes: tie down straps, orbit drive, all-terrain castors, flat stretcher top with tie down straps and gap bridge, fully customised top surface option to suit special purpose & comes standard without backrest

Options: Can be supplied with or without dropsides, emergency CPR backrest available

Product Code

Electric: MPT-C-HE

Electric Package: MPT-C-HE-P

Hydraulic: MPT-C-HH

v-17012024