

# Contour Classic

The most flexible custom design platform available for a transport and recovery stretcher.



## Key Benefits

**Customisable design:** Choose from a variety of dropside options and electric or hydraulic versions available.

**Effortless control:** Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

**Supercharged:** Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

**Custom quality:** Manufactured in Australia with speed, quality, and care, plus customisation capability.

**Longest warranty:** Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

## Application

Ideal for Day Surgery, Recovery, Transit Lounge & Cardiac

## Why Buy?

The most flexible platform of all our stretcher range. The twin column design allows for several height range choices - we can build in a low height range for patient egress care or higher height ranges for better equipment or surgeon access with less strain during procedures.

Tailor make the Contour Classic with shorter, longer, narrower and wider top sizes and select from a range of dropside styles to suit specialised procedures or the patient type you are working with. To simply put it, there is a high level of customisation possible with many accessory mounting options and solutions.

v-01032023

# Contour Classic



## Hydraulic option

Non-powered, manually operated stretcher option.

In some cases, such as small or remote facilities, or large high pressure departments, hydraulically operated stretchers may be more beneficial, with simple mechanical servicing and without the requirement of regular charging or power cords.



## Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 810mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



## Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.



## Ultra-low height

Very safe and easy for patients to get on and off.

The 440mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



## Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.

## Specifications

**Overall Size:** 2135 x 815mm (705mm narrow & 965mm wide options)

**Patient Surface:** 1900 x 680mm (570mm narrow & 830mm wide options)

**Electric Height Ranges:** 440 - 810 (std) up to 540 - 910mm

**Hydraulic Height Range:** 525 - 850mm

**SWL:** 300kg

**Castor Size:** 200mm castors and Ezi-Trak 5th wheel steering

**Electric Package Includes:** 100mm premia mattress, folding IV, vertical oxy, dropdown footboard, drop down push handles, 2nd handset & USB port

**Hydraulic Package Includes:** 100mm premia mattress, folding IV, vertical oxy, dropdown footboard & drop down push handles

## Product Code

Electric: MPT-C-CE

Electric Package: MPT-C-CE-P

Hydraulic: MPT-C-CH

Hydraulic Package: MPT-C-CH-P

v-01032023