

Contour Exami

An economical, customisable, easy to move examination table available with or without dropsides.



Key Benefits

Safe & low: The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting.

Excellent height: Electrically adjust the patient height up to 930mm helps to examine or treat your patient without having to stoop over, ensuring easy access and less stress.

Comfortable: 100mm fluid proof Protran mattress with a welded, pressure relieving cover.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for medical clinics, examination rooms, schools, prisons & mining sites

Why Buy?

Helping patients on and off an examination table is easy and safe with a design that allows the mattress platform down to 400mm, lower than a standard office chair height! Uniquely, these stretcher examination tables are designed to be easily moved around and re-positioned due to their light-weight build and large castors. There are a range of accessories and customisations available to ensure the unit is designed to be priced and configured best to practically fit into your modern facility.

v-10102023

Contour Exami



Without dropsides

If you require a patient examination table, then this is worth considering. The standard Contour Exami without dropsides gives all the same functionality as an examination table, however has the additional value of colour choices and accessory mounting, while being very easy to move and reposition as required.



Ultra-low height

Make it easier and safer for patients to get on and off. The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Excellent high height

Don't strain your back, raise a patient up to your best working height. Lifting the platform height up to 830mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release option available.



Head and foot tilt

A generous 15 degrees top tilt both directions. Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Customise it your way

Upgrades such as wider top, dropsides, base cover, colours and accessories. This is a very flexible platform with excellent customisation available. Change the standard Contour Exami from a basic examination table into an economy patient stretcher.



Go anywhere

Change from the standard 125mm twin to 175mm single castors for situations such as outdoor. If there is a requirement to wheel over rougher concrete, sand or grass terrains, this option gives larger diameter rubber wheels with 50mm extra ground clearance under the stretcher.

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 400 to 830mm

SWL: 300kg

Castor Size: 125mm twin wheel castors, 4 brakes and 1 steer

Includes: 100mm Protran mattress, electric backrest that raises to 85 degrees, standard main power electric operation

Options: dropsides, base cover, battery back-up with cable retractor, larger castors and 6 endrail colours

Product Code

Contour Exami (without dropsides); MPT-C-EXAMI

150mm wider top; MPT-CXA-WT

Dropsides; MPT-CXA-5RD

Base cover with storage; MPT-CXA-BC

Battery backup; MPT-CXA-BB

CPR backrest release; MPT-CXA-BQR

175mm large castor option; MPT-CXA-175

v-10102023